

Containment Swab™ (C-Swab)™

(Containment tube, longer handle to aid in collection of sexual assault evidence, avoid loss of sample material as the cotton tip is centrally positioned in the tube, reclosable, sterile, latex free, Made in USA)



Features and Benefits

- Tip:** High Quality USP grade cotton- carded and coiled onsite to maintain consistency. Bonded securely to a semi-flexible polystyrene handle. Diameter: 5.156 mm Length: 17.45 mm
- Handle:** White Semi-flexible polystyrene Tube. Diameter: 2.54mm. Length: 149.225 mm
- Tube:** Molded base holds tip in center of tube
- Uses:** **Prevents contamination, helps avoid loss of sample material as the swab tip is held in the center of the tube. Reclosable.**
- Standards:** ISO 13485:2003, ISO 9001:2008. Certificate of Sterility available with each shipment upon request.
- Qualities:** Absorbent tip. Transport device/swab as single unit. Sterile, Individually wrapped, Sealable by an evidence information label.
- Pack:** 50 per box, 400 per case (packed as 50/box i.e. 8 boxes) and 5000 per case (packed as 50/box i.e. 100 box)

containmentswab.com

www.c-swab.com

Abacus Diagnostics, USA
Dramatically Improving Crime Scene Investigation