


FAB-Swab Filtered Air Breathable Cotton Tipped Applicator



Features and Benefits

Tip: High Quality USP grade cotton- coiled onsite to maintain consistency. Diameter: 5.15mm/ 0.203". Length: 17.4mm/ .687"

Shaft: White polystyrene Breakable when bent to 90°. Diameter: 2.54mm/ 0.10". Length: 129mm/ 5.1"

Tube: Transport container with 8 built in micro-holes protected by a special breathable filter to prevent cross contamination during transport or storage. Length: 168mm/6.625"

Uses: Designed for convenient reliable sample and evidence collection. Can be used dry, with IPA, or sterile water.

Standards: ISO 13485:2003, ISO 9001:2000, FDA GMP's, Certificate of Sterility available with each shipment upon request.

Qualities: Sterile, Individually wrapped, Sealable by an evidence information label.

Pack: 50 per box, 400 per case (packed as 50/box i.e. 8 boxes) and 5000 per case (packed as 50/box i.e. 100 box)

Patent: Worldwide

www.abacusdiagnostics.com